



Astrocyte Thawing Medium

Catalog Number: MEDAST-T1

Volume: 50 mL

Product Format: Ready-to-use, sterile, serum-containing recovery medium

Lot Number: Provided on bottle label and Certificate of Analysis

PRODUCT DESCRIPTION

Astrocyte Thawing Medium (Catalog # MEDAST-T1) is a complete, sterile, optimized formulation specifically designed to maximize post-thaw viability, attachment, and functional recovery of cryopreserved human iPSC-derived astrocytes (e.g., Cat # IPSCDA001) and primary astrocyte cultures. The medium contains a proprietary combination of serum, antioxidants, and a cell-protective reagent that minimizes cryoprotectant (DMSO) toxicity, reduces osmotic shock, and prevents membrane damage during the thawing process.

Key Features:

- **High viability recovery:** ≥85% viable cells post-thaw (validated with Cat # IPSCDA001)
- **Serum-containing formulation** (10% FBS) – enhances membrane stabilization
- **Contains DNase I** – reduces clumping from dead cell DNA release
- **Contains antioxidant cocktail** – reduces post-thaw oxidative stress
- **Ready-to-use** – no additional supplementation required
- **Single-use format** – 50 mL sufficient for thawing up to 10 vials (5 mL/vial)

Note: Does **not** contain growth factors (FGF-2, EGF). For long-term culture, replace with Astrocyte Maintenance Medium (Cat # MEDAST-M1) 24 hours post-thaw.

STORAGE & HANDLING

Parameter	Condition
Storage temperature	2°C to 8°C (DO NOT FREEZE)
Shelf life (unopened)	3 months from date of manufacture
Shelf life (opened)	2 weeks when stored aseptically at 2–8°C
Protection from light	Not required (no light-sensitive components)
Pre-warming before use	37°C water bath for 10–15 minutes (DO NOT exceed 37°C)
Color indicator	Rosy pink (pH 7.2–7.4) → Yellow indicates pH drop (contamination or CO ₂ imbalance)
Aliquot recommendation	Do not refreeze; discard unused portion after 2 weeks

THAWING PROTOCOL (For Cat # IPSCDA001)

Preparation

1. **Pre-coat culture vessels** (ideally day before or at least 2 hours before thawing):
 - Smooth Coat Solution overnight at 37°C, rinse twice with PBS,
 - **Pre-warm** Astrocyte Thawing Medium to 37°C (water bath or incubator).

Thawing Procedure (Do Not Centrifuge Immediately)

- | Step | Action |
|-------------|--|
| 1 | Remove vial of astrocytes (Cat # IPSCDA001) from liquid nitrogen. |
| 2 | Loosen cap ¼ turn to relieve pressure, then retighten. |
| 3 | Thaw in 37°C water bath with gentle agitation until small ice crystal remains (~90–120 seconds). Do not over-thaw. |
| 4 | Decontaminate vial with 70% ethanol. |
| 5 | Immediately transfer cells dropwise into 10 mL pre-warmed Astrocyte Thawing Medium in a 15 mL conical tube. |
| 6 | Centrifuge at 300×g for 5–7 minutes at room temperature (not 4°C – cold reduces viability). |
| 7 | Aspirate supernatant (contains DMSO). |
| 8 | Gently resuspend pellet in 1–2 mL pre-warmed Astrocyte Thawing Medium . |
| 9 | Count cells (trypan blue exclusion). Expected viability: ≥85%. |
| 10 | Dilute to the desired seeding density in Astrocyte Thawing Medium (not maintenance medium yet). |
| 11 | Plate cells onto a coated vessel. |
| 12 | Critical step: 4–6 hours post-plating, replace Thawing Medium with Astrocyte Maintenance Medium (Cat # MEDAST-M1) to remove serum and DNase I. |

Why replace at 4–6 hours? Longer exposure to serum promotes astrocyte reactivity and fibroblast overgrowth. Short exposure (4–6h) is sufficient for attachment but minimizes phenotypic drift.



VALIDATION DATA (Using Cat # IPSCDA001)

Parameter	Result with MEDAST-T1	Result with standard DMEM/F12 + 10% FBS
Viability (immediate post-thaw)	88 ± 3%	72 ± 6%
Attachment efficiency (24h)	92 ± 4%	68 ± 8%
Clumping (visible aggregates)	Minimal (<5% of cells)	Moderate–high (15–25% of cells)
GFAP expression (day 5)	91 ± 4% positive	73 ± 9% positive
IL-6 response to TNF α /IL-1 β (day 7)	14-fold increase	6-fold increase

Data represent mean \pm SD from 3 independent lots of iPSC-astrocytes, n=6 replicates per condition.

QUALITY CONTROL SPECIFICATIONS



Test Parameter	Specification	Result (representative lot)
Appearance	Clear, pink solution	Clear, pink solution <input checked="" type="checkbox"/>
pH (at 25°C)	7.2 – 7.4	7.3 <input checked="" type="checkbox"/>
Osmolality (mOsm/kg)	320 – 360	338 <input checked="" type="checkbox"/>
Sterility (14-day culture)	No bacterial/fungal growth	Negative <input checked="" type="checkbox"/>
Mycoplasma (PCR)	Negative	Negative <input checked="" type="checkbox"/>
Endotoxin (LAL assay)	<1.0 EU/mL	0.4 EU/mL <input checked="" type="checkbox"/>
DNase I activity	Functional (degrades DNA gel)	Positive (30 min, 37°C) <input checked="" type="checkbox"/>
Cell performance test*	≥85% viability + ≥85% GFAP+ at day 5	88% viability, 91% GFAP+ <input checked="" type="checkbox"/>

*Performance tested using Human iPSC-Derived Astrocytes (Cat # IPSCDA001) thawed per protocol, plated on a Smooth Coat Solution Coated Plate.

COMPATIBLE CELL TYPES

Cell Type	Catalog #	Recommended?
Human iPSC-Derived Astrocytes	IPSCDA001	<input checked="" type="checkbox"/> Optimized
Human iPSC-Derived Astrocytes (P2)	IPSCDA002	<input checked="" type="checkbox"/> Yes
Primary human fetal astrocytes	Various	<input checked="" type="checkbox"/> Yes



Cell Type	Catalog #	Recommended?
Primary rat/mouse cortical astrocytes	N/A	☑ Yes
Human iPSC-Derived NPCs (for thawing)	IPSCN001	☒ No (use NPC Thawing Medium, Cat # MEDNPC-T1)
Human iPSC-Derived Neurons	IPSCNEU001	☒ No (use Neuron Thawing Medium)

TROUBLESHOOTING GUIDE

Problem	Possible Cause	Solution
Low viability (<70%)	Thawed too slowly; DMSO toxicity	Thaw rapidly (<2 min); use pre-warmed medium.
Excessive clumping	DNA release from dead cells	Medium contains DNase I; ensure 50 U/mL is active (store at 2-8°C).
Poor attachment (<50% after 24h)	Coating inadequate or cells seeded in maintenance medium	Use recommended coating; do not seed in serum-free medium.



Problem	Possible Cause	Solution
Floating cells after 4-6h	Normal (dead cells detach)	Gently wash with warm PBS before adding maintenance medium.
Cells look rounded at 24h	Serum reaction or endotoxin	Check endotoxin (<1 EU/mL); replace medium at 4-6h exactly.
Yellow medium before use	CO ₂ exposure during storage	Discard; do not use if pH is orange/yellow.

APPLICATIONS FOR RECOVERED ASTROCYTES

After thawing, astrocytes are suitable for:

Application	Time post-thaw
Immunocytochemistry (GFAP, S100 β , AQP4)	Day 5–7
Glutamate uptake assays	Day 5–7
Cytokine stimulation (TNF α , IL-1 β , LPS)	Day 6–8



Application

Time post-thaw

Neuron-astrocyte co-culture

Plate astrocytes day 0, add neurons day 5–7

Calcium imaging (spontaneous or ATP-evoked)

Day 6–10

Phagocytosis (myelin, synaptosomes)

Day 7–10

RNA extraction / qPCR

Day 5–10

PRECAUTIONS & LIMITATIONS

- **For research use only.** Not for clinical, diagnostic, or therapeutic use.
- **Do not use if medium is turbid, contains precipitate, or is yellow/orange.**
- **Do not freeze** – serum components will denature.
- **Do not use for neurons or NPCs** – serum in this medium induces neuronal differentiation and NPC adhesion changes.
- **Replace with maintenance medium within 4–6 hours** – prolonged serum exposure causes astrocyte reactivity and proliferation of non-astrocyte contaminants.
- **Use within 2 weeks of opening** – discard any unused medium after 14 days.

ORDERING INFORMATION & RELATED PRODUCTS



Product	Catalog Number	Size
Astrocyte Thawing Medium	MEDAST-T1	50 mL
Astrocyte Maintenance Medium (serum-free)	MEDAST-M1	500 mL
Astrocyte Maintenance Medium (low serum: 1% FBS)	MEDAST-M1-S	500 mL
Human iPSC-Derived Astrocytes (P1)	IPSCDA001	500,000 cells
Human iPSC-Derived Astrocytes (P2 expanded)	IPSCDA002	1 million cells
Coating Kit	ECM-AST1	1 kit
Coating Kit (for NPCs & astrocytes)	ECM-NPC1	1 kit
DNase I (animal-free, 2500 U)	ENZ-DN1	1 vial
Astrocyte Functional Assay Kit (IL-6 ELISA + Glutamate uptake)	ASY-AST1	96 tests

IMPORTANT: This warranty is valid for one month from the original purchase date.